

RNA Quality

Determine RNA quality with use of Agilent Bioanalyzer 2100.

Turbo DNA-free:

Add the following to a RNase-free, 0.5-ml microcentrifuge tube on ice:

Up to 10 μ g RNA

0.1 Volume Reaction Buffer

1 μ l Turbo DNase

Molecular Biology Grade Water to 10 μ l

Incubate 30 minutes at 37°C

Add 1 μ l Turbo DNase and incubate an additional 30 minutes at 37°C

Add 2 μ l of DNase Inactivation Reagent

Incubate 2 min. room temp.

Centrifuge 10,000 x g for 2 min.

Remove supernatant to new tube.

cDNA Synthesis Using SuperScript III

1. Prepare Master Mix

Oligo(dT) ₂₀	1 μ L
dNTP, 10 mM	1 μ L
RNA	1-10 μ g (no more than 10 μ L)
RNase-free Water	Up to 12 μ L

2. Heat mixture to 65°C for 5 minutes

3. Place on ice for 2 minutes.

4. Prepare Master Mix

5x First Strand Buffer	4 μ L
0.1 M DTT	1 μ L
RNaseOUT	1 μ L
SuperScript III	2 μ L

5. Mix by pipetting.

6. Incubate at 50°C for 60 min

7. Inactivate reaction at 70°C for 15 min

8. Cool to 4°C.

The cDNA can now be used as a template for amplification in PCR. However, amplification of some PCR targets (those >1 kb) may require the removal of RNA complementary to the cDNA. To remove RNA complementary to the cDNA, add 1 μ l (2 units) of *E. coli* RNase H and incubate at 37°C for 20 min.

Cost per reaction:

TurboDNase	\$3.32
Oligo(dT) ₂₀	\$2.00
dNTP, 10 mM	\$0.50
5x First Strand Buffer	-----
0.1 M DTT	-----
RNaseOUT	\$0.80
SuperScript III	\$9.00
RNase H	\$6.00

qRT-PCR Conditions and Analysis

PCR reactions are performed in an optical 384-well plate with an ABI PRISM 7900 HT sequence detection system (Applied Biosystems)

For 10 µl reaction volume:

2 µl of primer pair (1 uM)

2 µl of 1:20 dilution of cDNA

5 µl of 2X Power Sybr Green MASTER MIX (Cat no. 4368708)

1 µl Molecular Biology Grade Water

Turn on 7900HT

Open the “template 5uL” located in the microarray facility inbox

Change the Sample Volume to 10 µl

Click connect

Click Open/Close open door for placement of sample plate

Click Start

Data analysis:

Data is collected and analyzed using the SDS 2.2.1 software (Applied Biosystems).

Click the for automated data analysis.

Cost per reaction:

Power Sybr Green Master Mix	\$0.27
384-Well Plate	\$1.80
Optical Cover	\$1.00